

CardioChek® PA - Serial No:

Site Name:

Internal Quality Control Tests Log sheet - NB: For each CardioChek® PA device you must use a separate log sheet

Date	Reason for test - mark or tick (ONE only)			Control Solution & Level being tested				Lot number of control solution	Expiry date of control solution	Date of opening vial	Type of test strips	Batch number of test strips	Expiry date of test strips	Test	Expected results range (mmol/L)	Test result (mmol/L)	Pass (P)/ Fail (F)	Operator initials
	Routine	New strip batch	Suspected error	MC L1 (w)	MC L2 [r]	HDL L1 (y)	HDL L2 (g)											
01/01/2014	X	-	-	X	-	-	-	MC13A	28/02/2014	01/11/2012	TC/HDL Duo	U201	08/05/2013	CHOL	4.4 - 7.7	6.7	P	JAC